


1.3-inch Smart Heart Rate and Blood Pressure Bracelet

User Guide




Instructions for Quick Use of Bracelet

An Introduction to Components



Charge Mode

For the first time, please take out the bracelet and charge the device with the corresponding charger. (The charging method is as illustrated in the following figure: put the bracelet into the charger; aim it at two contacts of the charger, and cover the charger to charge. The charging voltage is 5V.)



Notifications

Press and hold the touch button for 3 seconds to check the historical information; If there's no historical information, it will be empty.



1.3-inch, large font, multi-display content, convenient check.



State


Switch to this interface to view the states of current steps, distances and calories. Distances and calories will be calculated and displayed according to the current step number, the height and weight set by the APP.

Notes for heart rate monitoring:
•The heart rate LED light behind the bracelet will be on when measuring.
•The sensor should be close to the skin to prevent external light from affecting the accuracy of the measurement;
Measuring area should be kept clean, because sweat or stains will affect the measurement results.



Way to wear

1. It is preferable to wear the bracelet behind the ulnar styloid process;
- 2.Adjust the size suitable for the wrist according to the regulation hole and buckle the wrist strap;
3. Keep the sensor close to the skin to avoid movement.



Switch on / off

- 1.When the bracelet is power off, charge it or press the function button for 3 seconds to start up the bracelet with shake;
2. When the bracelet is turned on, switch to the POWER OFF interface and press the function button for 3 seconds to shut down the bracelet with shake.

Bracelet Operation

1. When the bracelet is on, quickly touch the function key to light up the screen or switch the menu;
2. The screen will be off after 5 seconds without operations by default. The user can connect it to the phone APP to modify the screen off time;
- 3' The call and message reminder of the bracelet is open by default. The user can close it in the APP;
- 4.Switch it on or off in the main interface of (Android)/ Smart Play (IOS).

To install Bracelet APP

Scan the following QR code or enter the application markets to download and install it. "FlagFit2.0"



Device Requirements: IOS 8.0 and above versions; Android 4.4 and above versions; Support Bluetooth 5.0.

Device attachment


When being used for the first time, it is necessary to connect APP for alignment. After successful connection, wristband will self-synchronize time, or otherwise, step calculation and sleep data won't be accurate.

Open the phone APP. Click My Icon.

Add device (IOS) / Add device (Android)

Click the search icon in the upper right.

Click the device connection.



Bluetooth connection: after successfully matching, APP will automatically save the wristband bluetooth address. When APP is opened or operated in background, it will automatically search and link with the wristband.

Data Synchronism: manually pull down the synchronization data on the APP home page; the bracelet can store off-line data for 7 days, the more data, the longer the synchronization time is, and the longest time is about 2 minutes.

APP Functions and Settings

Personal Information and Sport Target

Please set the personal information after entering the APP.


My - Click the avatar - Avatar change/Background change/Cancel

My - Click the avatar - Avatar change/Background change/Cancel

- Personal settings-Sex-Age-Height-Weight, distance unit.
- Set personal daily sport target. Reasonable sports planning is conducive to the improvement of physical fitness.

Reminder Function

The reminder switch should be switched on the "Play Bar" page in the APP and keep the Bluetooth of the phone connected with the bracelet. It will shake if there's a message.



Call Reminder: When there is an incoming call, the bracelet shakes and displays the name or phone number of the caller. If the "Call Reject" is enabled on the APP, press and hold the function button for 3 seconds to reject the call.

To test your heart rate, you need to take a long time signal to accurately calculate your heart rate.

What do Bluetooth connections do? (connection / reconnection failure or connection speed is slow)

Confirm whether your device is 4.4 or more versions of IOS8.0 and Android system and support Bluetooth 4.

- Because of the problem of signal interference in Bluetooth wireless connection,
- a. there may be no time to connect each other. If the connection is not long enough, please ensure that the operation is carried out without magnetic field or without multi Bluetooth device interference.
- Turn off Bluetooth and open.
- b. Turn off the mobile phone background application or restart the mobile phone.
- c. Mobile phones do not connect to other Bluetooth devices or functions at the same time.
- APP is running normally in the background and may not be connected if it is not in the background.

Can't search the wristband?

Since the wristband is not searched by phone in Bluetooth broadcast, make sure the wristband is in power and activated and is not connected to any

Message Reminder: When receiving messages, the bracelet shakes and displays icon reminder. Click to view the message content.

QQ Reminder: When receiving QQ messages, the bracelet shakes and displays icons. Click to view the QQ content.

WeChat Reminder: When receiving WeChat messages, the bracelet shakes and displays icon reminder. Click to view the WeChat content.

Smart Alarm: A total of 3 alarms can be set, which, after being set, can be synchronized to the bracelet . Smart alarm is a silent alarm clock. When alarming, the bracelet shakes moderately and supports offline alarm reminder.

Sedentariness Reminder: The reminder interval is one hour by default. Under the connection state, if the user does not exercise within an hour, the bracelet will shake to remind the user to exercise.

Charge Reminder: When the bracelet is connected to the power supply or charging clip, the motor will shake. Please note that the bracelet will not shake when the power is low or full.

Sleep Detection

When you fall asleep, the bracelet will automatically judge and enter the

sleep monitoring mode, automatically detect the duration of your time of deep sleep/light sleep/wakefulness all night, and calculate your sleep quality. For the time being, sleep data can only be viewed on the APP. Attention: There is sleep data only when you go to sleep wearing the bracelet.

Main Functions

Main Interface

The bracelet has 3 main interfaces. Press and hold the main interface to switch the interface. After the APP is firstly matched and connected successfully, the bracelet will sync the phone's time, date and day of the week.



Does the Bluetooth connection need to be open all the time? Is there any data after disconnection?

Before the data is synchronized, the data is kept on the wristband host (seven days of data can be retained). When the Bluetooth connection of the wristband and the wristband is successful, the wristband data will be automatically uploaded to the mobile terminal. Please synchronize the data to the mobile terminal in time.

Note: if the caller and SMS alerts are opened, you must keep the Bluetooth connection, and the alarm is offline.

Basic parameter

| | | | |
|------------------|--------------------------------------|------------------|-------------------------------------|
| Type | Smart continuous heart rate bracelet | Battery type | Li-polymer |
| Display screen | 1.3 cun IPS | Battery model | 392121 |
| Weight | 31g | Battery capacity | 180mAh |
| Length | 25.0cm | Battery voltage | DC 3.7V |
| Waterproof level | Ip67 | Battery power | 0.666Wh |
| Material | Silica gel Watchband | Packing list | Host, instructions Charging clip |