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Android 4.2

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

# Smart Bracelet User Manual


**1. Compatible System of Phone and BLE Version :**

Android 4.2 and above

iOS 8.0 and above

BT4.0

## 2. User Preparation

Please make sure the device is full charged before the first time using.

Heart rate sensor  
USB interface

**3. APP "TFit" Down-load Scan the QR Code**

Scan the QR Code

**4. Bracelet Connection**

**5. Bracelet Operation**

**Smart Alarm:** Find the "Smart Clock" on APP, Set the time, press the time, can set the repeat mode, press complete to save the setting, the bracelet will vibrate when the time arrived.

**6. Function Instruction**

**7. Bracelet Operation**

**Heart rate monitor:** Automatically measure the heart rate when the function mode, press complete to save the setting, the bracelet will vibrate when the heart rate arrived.

**Pedometer:** 24 hours recording on bracelet

Clear bracelet data at 24:00, Can check data

Can press the heart rate monitor in APP, Data will be saved in the APP

**Sleep monitor:** Automatically record your sleep status with analyzing the deep sleep and light sleep hours, data saved in app.

**Sedentary Reminder:** Open the function in APP and set the time duration, the bracelet will vibrate

**Sports mode:** Statistics of your day Exercise duration, 24:00 Am per day View history in APP;

**Reminder for drinking:** Open the function in APP and set the time duration, the bracelet will vibrate

**Anti-lost:** Open the function in APP, the bracelet will vibrate when the phone is taken away

**SPO2 Test:** Automatically measure the SPO2 rate when the function switch to the Heart rate, or can press the heart rate monitor in APP, Data will be saved in the APP

**Message Reminder:** Open the function in APP, the bracelet will vibrate when the phone comes, a new message (SMS, QQ, Facebook, Twitter, WhatsApp, Wechat) will be saved in the APP

**Call Reminder:** Open the function in APP, the bracelet will vibrate when the phone comes

**Shake the bracelet to take photo:** Open the function in APP, Then press the "Shake to take a photo", enter into the camera, shake the bracelet to confirm. (Make sure the mobile phone setting: Allow the APP to take photo and save the photo)

**Find Bracelet:** Press the function in APP, the bracelet will vibrate.

**Restore factory Settings:** set the time to the previous day in the mobile APP Or the data of the bracelet can be cleared by synchronizing one day later;

\*The company reserves the right to make changes to the contents of the specification without notice.

APP Home Screen Examples:

1. Heart rate monitor screen: Shows Heart rate (86 BPM), Steps (08523), Calorie (0852), Mile (0.68), and Bpm (86).

2. Sleep monitor screen: Shows Sleep status (Sleep, Deep sleep, Shallow sleep) and Total time (00:00:00).

3. Sports mode screen: Shows Exercise duration (00:00:00) and Total time (00:00:00).

4. Anti-lost screen: Shows Blood oxygen (86 SPO2) and BP (120/80 mmHg).

5. Message reminder screen: Shows Blood oxygen (86 SPO2) and BP (120/80 mmHg).