









<h1 style="text-align: center;">Hello</h1> <h2 style="text-align: center;">User Manual</h2> <p style="text-align: center;">Thank you for choosing digital microscope</p>	<h2 style="text-align: center;">CONTENTS</h2> <p style="text-align: center;">English.....01-08</p> <p style="text-align: center;">Chinese.....09-16</p>	<h3>Preface</h3> <p>Thank you for choosing this equipment. The usage of this product is complicated, please read this manual carefully before using to understand the operation of microscope.</p> <h3>Important Notes</h3> <ol style="list-style-type: none"> <li>Please open the dust cover in front of the lens when using it for the first time.</li> <li>This device only supports Android mobile phones with OTG function. Please confirm whether the mobile phone supports OTG and check whether the OTG function of the mobile phone is turned on.</li> <li>The camera button on the product can only be used on the software "ViewPlayCap", other software cannot be used.</li> <li>This equipment is a combination of digital microscope and optical microscope. The specific magnification effect is subject to the actual picture taken.</li> <li>The best focal length of the microscope is 0-40mm. The image is not clear beyond 40mm, you need to adjust the focus by adjusting the focus wheel to see the object clearly.</li> </ol> <h3>Product introduction</h3> <p>Thank you for purchasing our WiFi electron microscope, which can be easily used in different areas including:</p> <ol style="list-style-type: none"> <li>Textile industry, used for textile inspection</li> <li>Printing inspection</li> <li>PCB, precision machinery testing</li> <li>Educational purpose</li> <li>Hair check</li> <li>Skin examination</li> <li>Microbial observation</li> <li>Inspection of jewelry and coins (collection)</li> <li>Visual AIDS</li> </ol> <p style="text-align: right;">- EN01 -</p>	<h3>Product description</h3>   <p style="text-align: center;">- EN02 -</p>	<h3>—Mobile users</h3> <p>Only supports Android phones with OTG function (Apple phones are not supported)</p> <h4>1. Software installation:</h4> <p>Open the browser and scan the following QR code to download and install the "USB CAMERA" software.</p>   <h4>2. Device connection:</h4> <p>Select the corresponding interface according to the type of mobile phone used, connect the device to the mobile phone, you can see the device LED will be on, at this time, if the device light is not on, please check whether the OTG function of your mobile phone is turned on or on the mobile phone settings page , or check if the interface is plugged in properly</p>  <p style="text-align: center;">- EN03 -</p>	<h3>—Computer users</h3> <h4>—Windows system:</h4> <p>Windows vista/7/8/10/11 or higher are supported</p> <h4>1. Software installation</h4> <p>Download and install the "ViewPlayCap" software, the download link is as follows: <a href="http://www.inskam.com/download/PC.zip">www.inskam.com/download/PC.zip</a></p>  <p>* When using the software "ViewPlayCap", the device camera button can be used.</p> <h4>2. Device connection:</h4> <p>a. Connect the USB cable of the device to the USB port of the computer, the light is on, and run "ViewPlayCap" b. Click on the device, switch and select the camera "HD camera",in the device, and you can use it, as shown in the figure below.</p>  <p style="text-align: center;">- EN04 -</p>	<p>c. To record video, you need to create a new folder, as shown in the following figure:</p> <p>Step 1: </p> <p>Step 2: </p> <p>Step 3: </p> <p style="text-align: center;">- EN05 -</p>	<h3>—MacOS system:</h3> <p>MacOS X 11.0 or later is supported.</p> <p>* When using the software "Photo Booth", the camera button of the device cannot be used, and the camera and video function can only be realized through the software</p> <ol style="list-style-type: none"> <li>Open the software (it is recommended to use the "Photo Booth" software that comes with the Apple computer).</li> <li>Connect the device</li> </ol> <p>1) First, in the "Applications" directory of the Finder window, find a software called "Photo Booth".</p>  <p>2. Connect the device</p> <ol style="list-style-type: none"> <li>Connect the device's USB cable to the computer's USB port, the light is on, and run "Photo Booth".</li> <li>Click on Photo Booth and select the camera "HD camera".</li> </ol>  <p>3. You can also use the alternative software "Quick Time Player".</p> <p style="text-align: center;">- EN06 -</p>	<h3>Frequently Asked Questions</h3> <p>Q1: What systems is this device compatible with? A1: Mobile phone support: Android 6.0+ (with OTG function) Computer support: Windows Vista/7/8/10/11+, Mac OS X 11.0+ Linux/Chrome OS is not supported.</p> <p>Q2: I just got it and connected it to my computer, it looks like the camera button doesn't work? A2: The camera button on the product can only be used on the software "ViewPlayCap", other software cannot be used.</p> <p>Q3: Why is the image very blurry when viewing? A3: 1. Have you adjusted the microscope to the best observation distance? You can adjust the focus by slowly rotating the focusing wheel left and right to achieve the clearest state. 2. Whether to open the dust cover in front of the lens.</p> <p>Q4: What is the reason for the unmoving black patches displayed on the screen when using the device to watch? A4: Due to some reasons, dust has entered. You can try to tap the product to see if the dust can be knocked out. If not, return it to the factory for repair.</p> <p>Q5: The microscope is connected to the PC, after successfully downloading and opening the software, there is no image. A5: 1. Reconnect the microscope and PC side. 2. Connect to the USB port on the back of the computer host.</p> <p style="text-align: right;">- EN07 -</p>	<p>Q6: After the device is successfully connected, what is the reason why the LED light does not light up and the display has no image? A6: 1. Confirm whether the light is turned off. 2. Reconnect the microscope to the computer or mobile phone. 3. Check whether the OTG function is turned on in the phone settings.</p> <p>Q7: After the device is successfully connected, the LED light is on and the phone has no image. What is the reason? A7: The mobile phone has OTG function, but without UVC image protocol, the microscope does not support this mobile phone.</p> <p>Q8: Can the device record sound? A8: No, it cannot record sound.</p> <h3>Technical Parameters</h3> <table border="1"> <tr><td>Pixel</td><td>0.3/2 million pixels</td></tr> <tr><td>Magnification</td><td>500X/1000X/1600X</td></tr> <tr><td>Picture resolution</td><td>640*480P/1280*720P/1920*1080P</td></tr> <tr><td>Focal distance</td><td>Manual focus(0~40mm)</td></tr> <tr><td>Picture format</td><td>JPG</td></tr> <tr><td>Video format</td><td>MP4/AVI</td></tr> <tr><td>Dynamic frame rate</td><td>20-30f/s</td></tr> <tr><td>USB interface type</td><td>USB 3.0/2.0/1.0</td></tr> <tr><td>Supportive system</td><td>Android 6.0+ Windows vista/7/8/10/11+ , Mac OS X 11.0+</td></tr> </table> <p style="text-align: right;">- EN08 -</p>	Pixel	0.3/2 million pixels	Magnification	500X/1000X/1600X	Picture resolution	640*480P/1280*720P/1920*1080P	Focal distance	Manual focus(0~40mm)	Picture format	JPG	Video format	MP4/AVI	Dynamic frame rate	20-30f/s	USB interface type	USB 3.0/2.0/1.0	Supportive system	Android 6.0+ Windows vista/7/8/10/11+ , Mac OS X 11.0+
Pixel	0.3/2 million pixels																										
Magnification	500X/1000X/1600X																										
Picture resolution	640*480P/1280*720P/1920*1080P																										
Focal distance	Manual focus(0~40mm)																										
Picture format	JPG																										
Video format	MP4/AVI																										
Dynamic frame rate	20-30f/s																										
USB interface type	USB 3.0/2.0/1.0																										
Supportive system	Android 6.0+ Windows vista/7/8/10/11+ , Mac OS X 11.0+																										

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