



Test Report

No. BOYRUWYR849217L1a

Sample Description

Disposable medical Mask

Applicant

Chongqing COE Display Technology Co., Ltd

Pony Testing International Group
www.ponytest.com



Statement

- 1.This report is invalid without the approver's signatures and special seal of inspection.
- 2.The pattern and characters of "PONY" and "谱尼" used in this report are protected by the trademark law of the People's Republic of China.Any unauthorized usage, counterfeit, forgery and alteration of trademarks of "PONY" and "谱尼" are the violations of the law. The PONY has the right to pursue all legal liabilities of the subject of the delict.
- 3.If the applicant has any questions about the results, shall provide a written retest application, the original report and prepay the retest fees to PONY within fifteen days since the approval date.
- 4.After the applicant finishes the procedure mentioned above, PONY shall arrange the retest as soon as possible. If the retest result accords with the applicant dissent, PONY shall refund the retest fees.
- 5.Tests that can not be repeated and tested shall not be carried out again.
- 6.The applicant should undertake the responsibility for the provided samples' representativeness and document authenticity. Otherwise, PONY has not any relevant responsibilities.
- 7.This report is only responsible for the test results of the tested sample. The test data only represents the evaluation of the tested sample and can only be used for researching, teaching or internal quality control by the applicant. PONY shall not be responsible for any economical or legal liability generated from direct or indirect usage of the test result.
- 8.PONY has the right to dispose the tested sample after approval of the test report.
- 9.PONY assures objectivity and impartiality of the test, and fulfills the obligation of confidentiality for applicant's commercial Information, and technique document.
- 10.The report is invalid in case of illegal transfer, embezzlement, imposture, modification or any altering, reproducing except in full, without approval of PONY. PONY shall investigate and affix the applicant's legal liability accordingly.

▲ Anti-counterfeiting Description:

- (1) The test report has exclusive report code.
 - (2) The test report is printed by anti-copying paper whose surface shows "PONY" security print with specific anticounterfeiting technique. Security print will disappear after copying. Duplicates are not expected to give "PONY" security print under any circumstances.
-



Hotline
400-819-5688

WWW.PONYTEST.COM

Scan QR code to follow
Wechat public account
"PONY4008195688"




Beijing Lab: (010) 83055000	Wuhan Lab: (027) 83997127	Harbin Lab: (0451) 58627755
Shanghai Lab: (021) 64851999	Changchun Lab: (0431) 85150908	Shijiazhuang Lab: (0311) 85376660
Qingdao Lab: (0532) 88706866	Dalian Lab: (0411) 87336618	Urumchi Lab: (0991) 6684186
Shenzhen Lab: (0755) 26050909	Zhengzhou Lab: (0371) 69350670	Hohhot Lab: (0471) 3450025
Tianjin Lab: (022) 23607888	Xi'an Lab: (029) 89608785	Hangzhou Lab: (0571) 87219096
Suzhou Lab: (0512) 62997900	Taiyuan Lab: (0351) 7555762	Ningbo Lab: (0574) 87736499
		Wenzhou Lab: (0577) 88271060
		Hefei Lab: (0551) 63843474
		Guangzhou Lab: (020) 89224310
		Xiamen Lab: (0592) 5568048
		Chengdu Lab: (028) 87702708

Test Results

No. BOYRUWYR849217L1a

Page 1 of 2

Sample Description	Disposable medical Mask	Sample Specification	Specification: 17.5*9.5 model: CCMM-SA sterile
Applicant	Chongqing COE Display Technology Co., Ltd	Trade Mark	---
Received Date	2020-04-02	Manufacturing Date or Lot No.	Manufacturing Date: 2020.3.17 Manufacturing Lot No.: 200317-C-S-01
Test Date	2020-04-02~2020-04-17	Sample Grade	---
Sample Quantity	50pcs	Test Type	Commissioning Test
Sample Status	Bulk (simple packaging)	Test Environment	To meet the requirements
Test Items	See the next page		
Reference Methods	See the next page		
Main Instruments	Electric Heating Constant Temperature Incubator, Respirator particle filtration efficiency and airflow resistance tester etc.		
Note	Limit on: EN 14683: 2019 This test report is to replace the test report No.BOYRUWYR849217L1. The No. BOYRUWYR849217L1 test report is invalid and of no legal effect. All related information should be referred to the new test report.2020-04-18 The information of above sample was provided by the applicant.		
	Edited by	Bian Youjia	
	Checked by	Guo Jun	
	Approved by	Da: GJ	
	Issued Date	2020-04-18	

Test Results

No. BOYRUWYR849217L1a

Page 2 of 2

Sample Description and Number	Test Items	Limit	Test Results		Evaluation	Reference Methods
		Type I				
R849217L1 Disposable medical Mask	Bacterial filtration efficiency(BFE), %	≥ 95	No.1	98.9	Conform	EN 14683: 2019
			No.2	97.6		
			No.3	98.9		
			No.4	98.0		
			No.5	96.9		
	Differential pressure, Pa/cm ²	<40	28.0		Conform	EN 14683: 2019
Microbial cleanliness, cfu/g	≤ 30	<1		Conform	EN 14683: 2019	

Sample Number and Photo:



Pony authenticate the photo on original report only

—End of Report—